


**Mykroscope**  
LEARNING TOGETHER

# Mykroscope Build Component

- 1 Laser cut cardboard ( 1.5mm )



- 2 Laser cut transparent acrylic ( 2mm )



- 3 Rocket switch ( 15 \* 10 mm )

- 4 AA battery holder ( 2 x AA Batteries )

- 5 LED light panel ( 37 \* 26 \* 2.3 mm )

- 6 Magnifying lens

- 7 AA batteries x 2

- 8 JST connector female & male ( optional )



**Mykroscope**

Emerging as champion of the Young Innovators Challenge 2017 (YIC) has enabled Sekolah Agama Menengah Jeram to turn their idea into reality. Their idea, the Mykroscope, is the pairing of a portable microscope with a smartphone. Mykroscope allows students to capture their discoveries on their mobile phones and share them with others around the world. With Mykroscope, students are no longer bound by the confines of their classrooms; the world is now their classroom.

[www.mykroscope.com](http://www.mykroscope.com)

#### Inventors

- Nor Asikin Bt Zulkifli
- Ainurania Azzahra Bt Zulkifli
- Nut Anis Bt Suratiman

#### Teacher Advisor

- Pn. Sabariah Binti Hussein

#### Mentor

- Muhammad Amir Zahidin Bin Zainal Abidin


#### Student of



SAM Jeram


Managed by  
**chmbaka**

# Mykroscope Manual





1



Oriентate and place the Mykroscope as shown on a flat table.




2

Fold the Mykroscope.









**8** Fold the Mykroscope.

**9** Insert the Lens.




**10** Fold the Mykroscope.


**11** Fold the Mykroscope.




**1 1** Fold the Mykroscope.




**1 2** Fold the Mykroscope.



**1 3** Fold the Mykroscope.




**1 4** Fold the Mykroscope.



1 5

Zip the Cable tie around the Slide Holder handle (Both side).



1 6

Push the Cable Tie close to the cardboard.